

Westfield Medical Ltd

Nominal Dimensions of Self-Seal Pouches

Product Code	Quoted Size	Actual Sizes (+/-3mm)					
		A	B	C	D	E	F
WPSS060100	60 x 100	57	127	100	20	2	7
WPSS090140	90 x 140	85	168	140	20	4	7
WPSS090200	90 x 200	85	230	200	20	3	11
WPSS090230	90 x 230	85	257	230	20	6	7
WPSS135255	135 x 255	134	282	255	20	5	7
WPSS135360	135 x 360	134	367	340	20	5	9
WPSS190330	190 x 330	190	360	330	20	5	9
WPSS305381	305 x 381	380	333	305	20	2	9

Values shown in shaded cells are critical dimensions, all other dimensions are shown for information only and do not form a specification.

Heat Seal Width = 10mm nominal, 6mm minimum.

Tape Width = 15mm

Print Position = 2mm (+3/-1mm) from end of sealing flap.

Note: Alternate chevron design may be used.

Westfield Medical Ltd
Nominal Dimensions of Self-Seal Pouches

